

PrestoCraft 3D

SHAPING

TOMORROW



DENTAL CASTABLE RESIN

Dental Castable Resin, a technological marvel meticulously designed to elevate the accuracy and intricacy of dental restorations. Tailored for high-resolution 3D printing, this resin redefines the casting process, allowing dental professionals to achieve unparalleled precision in the creation of intricate patterns for crowns, bridges, and a myriad of dental prosthetics.

KEY FEATURES:

- Negligible Ash Content
- Good Toughness
- Easy BurnOut
- Clean Precise Casting Results

Flexural Strength (MPa)	60-80
Flexural Modulus (MPa)	1500-2000
Residual Ash Content	<0.1%



Scratch
Proof



Light
Odor



Fast
Curing



Less
Brittle



Higher
Tenacity




Consistent
Quality

DENTAL RESIN

PrestoCraft 3D

DENTAL RESIN



DENTAL MODEL + RESIN

Step into the future of dentistry with our revolutionary Dental 3D Printing Model Resin, meticulously crafted to elevate precision, efficiency, and innovation in dental applications. This advanced resin is designed to meet the exacting demands of dental professionals, offering a seamless blend of cutting-edge technology and unmatched material quality.

KEY FEATURES:

- Fabrication of 3D Crown and Bridge Models
- Orthodontic Models
- Diagnostic Models
- Smooth Surface with Hard Finish
- Best Choice for Making Aligners

Biege

Grey



Consistent
Quality



Higher
Tenacity



Less
Brittle



Fast
Curing



Light
Odor



Scratch
Proof

DENTAL MODEL RESIN

PrestoCraft 3D

Achieve precision and excellence in Dental Modelling using our cutting-edge resin formulations. Crafted to produce detailed and accurate dental models for a variety of applications.

KEY FEATURES:


- Exceptional Detail and Precision
- Consistent Performance

PERFECT FOR:

- Orthodontic Study Models
- Prosthodontic Models
- Getting Perfect Fits for lab analogue
- Aligners making

Biege

Grey



DENTAL RESIN



Scratch
Proof



Light
Odor



Fast
Curing



Less
Brittle



Higher
Tenacity



Consistent
Quality

SURGICAL GUIDE RESIN

DENTAL RESIN

A dental 3D printing surgical guide is a specialized tool used in dentistry to precisely place dental implants with accuracy and efficiency.

Guide Placement: During the implant surgery, the surgical guide is placed over the patient's teeth or jawbone.

KEY FEATURES:

- Accurate Placement
- Streamlined Workflow
- Minimized Discomfort
- Predictable Results
- Tailored Design
- Avoids Errors

Colour	Clear-translucent
Wavelength	385-405
Odour	Low
Packaging	500g Bottle Net



Accurate



Biocompatible



Low Viscosity



Fast Curing

C&B CERAMIC RESIN

Prestocraft C&B ceramic resin offers biocompatibility, exceptional aesthetics, and robust strength, making it ideal for long-term prosthetic applications.

KEY FEATURES:

- Exceptional Strength & wear Resistance
- Exceptionally high aesthetic results after polishing
- Biocompatible Tooth colour Resin for long term Prostheses

Shades	BL, B1, A1, A2, A3
Wavelength	385-405
Odour	Low
Compressive strength	> 225mpa
Flexural Strength	> 150mpa
Packaging	500 g Bottle Net

Strength values depend upon adhering to recommended printing and post curing protocol

PrestoCraft 3D



Accurate



Biocompatible



Low Viscosity



Fast Curing

PrestoCraft 3D

CROWN & BRIDGE RESIN

Prestocraft C&B resin is a biocompatible, highly aesthetic formulation with inert fillers, designed to mimic the natural look of teeth.

KEY FEATURES:

- High Strength
- High Polishability for Tooth like Aesthetics
- Biocompatible Tooth Coloured filled Resin for Provisionals

Shades	BL, B1, A1, A2, A3
Wavelength	385-405
Odour	Low
Compressive strength	> 180mpa
Flexural Strength	> 100mpa

Strength values depend upon adhering to recommended printing and post curing protocol



DENTAL RESIN



Accurate




Biocompatible



Low Viscosity



Fast Curing



PrestoCraft 

Orthoplast

Near Ram Mandir,
Junction Road Khurja -203131
Uttar Pradesh (INDIA)
Phone: +91 94128 30836
GST No : 09AACFO2767J1ZX