



Premium Products for Premium Dentures



Scan the QR code
for more information



Prestorock Upper Anterior

Solvent Resistant Cross Linked
Double Acrylic Teeth

Prestorock Acrylic teeth are produced from Poly Methyl methacrylate in two chromatic layers

Shades : A1, A2, A3, A3.5, B1, B2

Moulds: S13, S14, S15, S18, S25, S32, S33, S38, S39, S41, S42, S51, S56, S67

Packing : A box contain 16 strips

Prestorock Lower Anterior


Solvent Resistant Cross Linked
Double Acrylic Teeth

Prestorock Acrylic teeth are produced from Poly Methyl methacrylate in two chromatic layers.

Shades : A1, A2, A3, A3.5, B1, B2

Moulds : D 35, D 37, D 39, D 45

Packing : A box contains 16 stripes



Prestorock Upper & Lower Posterior


Solvent Resistant Cross
Linked Double Acrylic Teeth

Prestorock Acrylic teeth are produced from Poly Methyl methacrylate in two chromatic layers.

Shades : A1, A2, A3, A3.5, B1, B2

Moulds : D 35, D 37, D 39, D 45.

Packing : A box contains 10 stripes (5 upper + 5 lower)





Prestorock Acrylic teeth Full Set

Solvent Resistant Cross Linked Double Acrylic Teeth

Prestorock Acrylic teeth are produced from Poly Methyl methacrylate in two chromatic layers.

Shades : A1, A2, A3, A3.5, B1, B2 **Moulds:** S13, S14, S15, S18, S25, S32, S33, S38, S39, S41, S42, S51, S56, S67

Packing : A box contains 4 Full Sets



Acryton Modelling Wax

Acryton dental modelling wax is for making occlusion rims and for setting up artificial teeth like PRESTOROCK during the fabrication of dentures.

Color : Pink | **Packing :** 12 sheets

Acryton Cold Mould Seal (Seprating Medium)

Aqueous solution can make strong seprating film on the surface of plaster model. Cold mould seal is a specially developed formulation of seprating solution for universal use, suitable for both packing/pressing technique and cast moulding.

Shades : Pink.

Packing : 110 ml Plastic Bottle, 450ml Plastic



Acryton Cold Cure Powder ►

An adaptable and aesthetically pleasing denture is achievable without any air bubble for motion or shrinkage. Variety of shade, Ultra fine particles promote superior burnishing & polishing.

Shades : Special Pink, Pink, Light Pink, Clear
Translucent, Light Veined, Veined

Packing : 110gm, 400gm, 3kg, 30kg.

Lab Pack : 400gm. Powder + 400ml. Liquid.



Prestocryl Heat Cure Powder with Hi-Strength ►

An adaptable and aesthetically pleasing denture is achievable without any air bubble for motion or shrinkage, Variety of shade, Ultra-fine particles with high strength promote superior burnishing & polishing.

Shades : Light veined, veined

Packing : 450 gms, 3 kg, 30 kg



Acryton Rapid Repair Liquid ▶

Acryton Rapid Repair Liquid binds with powder to help create a perfect denture with perfect shade and great hardness.

Packing : 110ml, 400ml, 4ltr, 40ltr, 200ltrs.



◀ Acryton Heat Cure Liquid

Acryton Heat Cure Liquid binds with powder to help create a perfect denture with perfect shade and great hardness.

Packing : 110ml, 400ml, 4ltr, 40ltr, 200ltrs.



Acryton Soft Liner ▶

A Temporary Soft - Liner is placed in a new or old denture in order to help improve the health of your gum tissues by absorbing some of the pressures of mastication.

- Useful for loose & irritating denture.
- Increase comfort during chewing and daily wear.
- Saves your gums from wear & tear of your hard plastic denture.

Shades : White

Packing : 56gms Powder and 56ml Liquid



Acryton Self Cure Temporary Crown Material



Acryton Self Cure Temporary Crown Material delivers superior handling properties, increased solvent resistance and higher melt temperature. It also has a significantly higher molecular weight range than standard tooth polymers.

Shades : A1, A2, A3

Packing : 56 gms , 400 gms

Acryton High Impact Heat Cure Powder

Acryton High Impact Heat Cure Powder provides additional strength to the denture from falls and other impacts while maintaining the perfect shades of denture.

Shades : Pink

Packing : 450gms



CERTIFICATE OF COMPLIANCE



INTERNATIONAL CERTIFICATION SERVICES PVT. LTD.
This is to certify that the
QUALITY MANAGEMENT SYSTEM of
ORTHOPLAST
Junction Road, Near Ram Mandir, Khurja 203131 Uttar Pradesh, India.
has been assessed and registered as complying with the requirements of the following International Standard:

ISO 13485 : 2016
The Quality Management System for Medical Devices applicable to :

Scope : Manufacturing and Supply of Artificial Teeth Sets, Denture Base Materials, Cold Mould Seal, Restorative Materials, ETCHANT, BOND, Temporary Crown and bridge material, Denture Relining Material.

Accredited by National Accreditation Board For Certification Bodies, India.

Validity of this certificate is based on periodic audits of the management system defined by the above scope and is contingent upon prompt, written notification of significant changes to the management system and/or its components thereof shall be immediately communicated to ICS. Further clarifications regarding the scope of this certificate and the applicability of ISO 13485:2016 requirements may be obtained at www.ics98888.com or www.ics98888.com

22/23, Goodwill Premises, Swastik Estate, 178 CST Road, Kallina, Santacruz (E), Mumbai - 400 098, Maharashtra, India. Tel: 022-42200800



Government of India
Ministry of Commerce and Industry
Directorate General of Foreign Trade
Office of the Joint Director General of Foreign Trade, Kanpur
4/286, Parvati Bagia Road, Kanpur

Importer-Exporter Code

This is to certify that **ORTHOPLAST** is issued an Importer-Exporter Code (IEC) **0513055045** with details as follows -

IEC	0513055045
रजिस्ट्रार का नाम (Firm Name)	AACFO2767J
फर्म का नाम (Firm Name)	ORTHOPLAST
विषय की प्रकृति (Nature of Concern)	Partnership
नियुक्ति की तिथि (Date of Issue)	25/10/2013
पता (Registered Address)	NEAR RAM MANDIR, JUNCTION ROAD Contact No: 9412641804, KHURJA UTTAR PRADESH, BULANDSHAHAR, UTTAR PRADESH, 203131
आंक का नाम / Name of the Signatory	GOVIND TAYAL
Director / Partner Details	Refer online at https://dgft.gov.in or scan the QR Code
शाखा/दफ्तर / Branch Details	Refer online at https://dgft.gov.in or scan the QR Code

Last Modified : 29/03/2022
File Number : KNPIECAMEND00012409AM22



Note : This is a system-generated certificate. Authenticity / Updated details of the IEC can be checked at official DGFT website <https://dgft.gov.in> by entering the IEC and Firm Name under Services > View Any IEC Details. You can also authenticate the certificate by scanning the QR code.

प्रेषक,
औषधि अनुज्ञापन एवं निर्वहन प्राधिकारी,
कार्यालय आयुक्त,
खाद्य सुरक्षा एवं औषधि प्रशासन, उ०प्र०,
सेक्टर-सी, अलीगंज, लखनऊ।

सेवा में,
M/s Orthoplast
Near Ram Mandir, Junction Road, Khurja,
District- Bulandshahar, Uttar Pradesh (India) - 203131

सं: MFG/MD/2023/83787 लखनऊ दिनांक 13.04.2023

विषय : प्रपत्र-MD-5 पर मेडिकल डिवाइस निर्माण अनुमति प्रदान किये जाने के सम्बन्ध में।

महोदय,
उपयुक्त विषयक आपके ऑन-लाईन आवेदन पत्र पर विचारोपरान्त कार्य एमडी-5 पर मेडिकल डिवाइस निर्माण हेतु लाइसेंस संख्या- MFG/MD/2023/000377, दिनांक 13.04.2023 मेडिकल डिवाइस रूल्स, 2017 में निहित प्रावधानों एवं औषधि एवं प्रसादन सामग्री अधिनियम, 1940 के अमीन स्वीकृत कर प्रेषित किया जा रहा है।

संलग्नक-उपरोक्तानुसार।

भवदीय

(सहायक चीफ़) **औषधि अनुज्ञापन एवं निर्वहन प्राधिकारी, उत्तर प्रदेश।**

सं: MFG/MD/2023/83787 तद्विनांक
प्रतिनिधि निम्नलिखित को सूचना एवं आवश्यक कार्यवाही हेतु प्रेषित-
(1) सहायक आयुक्त (औषधि), खाद्य सुरक्षा एवं औषधि प्रशासन, कार्यालय मण्डलानुसार, मेरठ मण्डल।
(2) औषधि निरीक्षक, खाद्य सुरक्षा एवं औषधि प्रशासन, कार्यालय जिलाधिकारी, बुलन्दशहर।

SANJIVA KUMAR CHAURASIA Digitally signed by SANJIVA KUMAR CHAURASIA
Date: 2023.04.13 14:02:44 +05'30'

(सहायक चीफ़) **औषधि अनुज्ञापन एवं निर्वहन प्राधिकारी, उत्तर प्रदेश।**



Please scan the above QR code to
visit our website and know more



ORTHOPLAST
NEW INDIA DENTAL PRODUCTS

Junction Road Industrial Area, Near Ram Mandir, Khurja (U.P.)-203131, INDIA

Email : nidproducts@yahoo.com | Website : www.nidproducts.com

Customer Care No.: +91-9412641804, 9219177654